

SoftGrip™ Pipette Kits

The best performance and comfort at a great value!

Hamilton's SoftGrip™ Manual Pipettes are easy to use, reliable, and maximize the efficiency of your lab. With three kits to choose from, you can get the best pipetting performance and comfort at an even better value.



Adjustable Volume Pipette Kit

- Any three SoftGrip™ adjustable volume pipettes
- Receive three racks of matching AdvanTip™ pipette tips
- 8-Position pipette carousel
- Free shipping

\$599

Value up to \$1,100

Fixed Volume Pipette Kit

- Any three SoftGrip fixed volume pipettes
- Receive three racks of matching AdvanTip pipette tips
- 8-Position pipette carousel
- Free shipping

\$399

Value up to \$790



Multi-Channel Pipette Kit

- Any SoftGrip 8-channel or 12-channel pipette
- Any box (10 racks) of AdvanTip pipette tips
- One pipette stand
- Free shipping

\$599

Value up to \$940

Call **1-888-525-2123** to learn more
or visit www.hamiltoncompany.com/pipettekits

Build Your Custom Kit

Use the tables below to select the pipettes in your kit. For more information on Hamilton pipettes, visit www.hamiltoncompany.com/pipettes.

Visit www.hamiltoncompany.com/pipettekits or call 1-888-525-2123 to learn more about SoftGrip pipette kits or to place an order.

Single Channel Adjustable Volume Pipettes

P/N	Volume Range	Color	Increments (µL)	At 10% of Pipette Volume		At 100% of Pipette Volume	
				Accuracy within ± %	Precision within + %	Accuracy within ± %	Precision within + %
55019-30	0.2–2 µL	Aqua	0.002	8.00	4.00	1.20	0.60
55019-32	1–10 µL	Purple	0.02	2.50	1.20	0.80	0.40
55019-34	2.5–25 µL	Forest Green	0.02	4.50	1.50	0.80	0.20
55019-36	10–100 µL	Violet	0.2	1.80	0.70	0.60	0.15
55019-38	30–300 µL	Brick Red	0.2	1.20	0.40	0.40	0.15
55019-40	0.100–1 mL	Sky Blue	2.0	1.60	0.50	0.30	0.12

8-Channel Adjustable Volume Pipettes

P/N	Volume Range	Color	Increments (µL)	At 10% of Pipette Volume		At 100% of Pipette Volume	
				Accuracy within ± %	Precision within + %	Accuracy within ± %	Precision within + %
55019-50	5–50 µL x 8	Sandstone	0.02	2.00	2.00	0.80	0.35
55019-52	30–300 µL x 8	Brick Red	0.2	1.20	0.70	0.50	0.30

12-Channel Adjustable Volume Pipettes

P/N	Volume Range	Color	Increments (µL)	At 10% of Pipette Volume		At 100% of Pipette Volume	
				Accuracy within ± %	Precision within + %	Accuracy within ± %	Precision within + %
55019-54	5–50 µL x 12	Sandstone	0.02	5.00	2.00	1.00	0.50
55019-56	30–300 µL x 12	Brick Red	0.2	2.40	1.20	1.00	0.50

Single Channel Fixed Volume Pipettes

P/N	Volume Range	Color	At 100% of Pipette Volume	
			Accuracy within ± %	Precision within + %
55019-01	5 µL	Steel Blue	1.40	0.75
55019-03	10 µL	Purple	0.80	0.40
55019-05	25 µL	Forest Green	0.80	0.30
55019-07	50 µL	Sandstone	0.60	0.30
55019-09	100 µL	Violet	0.50	0.20
55019-11	200 µL	Mustard	0.40	0.18
55019-13	250 µL	Burnt Orange	0.40	0.18
55019-15	300 µL	Brick Red	0.40	0.18
55019-17	500 µL	Olive Green	0.40	0.18
55019-19	1 mL	Sky Blue	0.30	0.12

Accessories

P/N	Product	List Price
6600849-01	8-Position Pipette Carousel	\$140.00
34195	Multi-Channel Pipette Stand	\$77.30



©2018 Hamilton Company. All rights reserved.
All other trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.
Lit. No. L70068 Rev D — 10/2018

HAMILTON®

Web: www.hamiltoncompany.com
USA: 800-648-5950
Europe: +40-356-635-055

Hamilton Americas & Pacific Rim
Hamilton Company Inc.
4970 Energy Way
Reno, Nevada 89502 USA
Tel: +1-775-858-3000
Fax: +1-775-856-7259
sales@hamiltoncompany.com

Hamilton Europe, Asia & Africa
Hamilton Central Europe S.R.L.
str. Hamilton no. 2-4
307210 Giarmata, Romania
Tel: +40-356-635-055
Fax: +40-356-635-060
contact.lab.ro@hamilton-ce.com

To find a representative in your area, please visit hamiltoncompany.com/contacts.