



HAMILTON COMPANY
P.O. BOX 10030
RENO, NEVADA 89520-0012 U.S.A.
TELEPHONE +1-775-858-3000
TOLL-FREE 800-648-5950
FAX +1-775-856-7259

Trueness and Precision Statement of Conformance

Hamilton syringes are manufactured to have Trueness (average accuracy) within +/- 1% of nominal volume and Precision (Coefficient of Variation) within 1%, measured at 80% of total scale volume. Design of our barrel and plunger assures high levels of Accuracy and Precision. The Hamilton Company Quality System is certified to ISO 9001 and FDA regulated. As such, our manufacturing processes, inspections and test methods utilize equipment and instrumentation that are calibrated through an unbroken chain of calibration, each traceable to NIST. The Hamilton Company calibration system conforms to the requirements of ANSI Z540.3

A handwritten signature in black ink, appearing to read "Eric Lute".

Eric Lute
Director of Quality

10-Feb-2017